

Beating Freeness Tester



Product Features

- The instrument has an adjustment mechanism
- The instrument uses manual compression to seal the vertebral body
- Use a three-legged leveling mechanism to adjust the centering of the sealed vertebral body to ensure the tightness of the sealed vertebral body and cylinder

Executive Standard

ISO 5267-1

Technical Parameters

No.	Parameter item	Technical index
1	Measuring range	0 ~ 100 S
	e	R
2	Ascending speed of the sealing cone	(100 ±10) mm /s
3	Bottom hole flow rate	(149±1) s
4	Remaining volume	7.5 ~ 8.0 mL
5	Dimension	450×250×1100 mm
6	Net Weight	30kg around